Solution for Genetic Technology

# **SolGent** Molecular Diagnostics

Real-Time RT-PCR for Detection **2019-Novel Coronavirus**

*Approved at third by KCDC-EUA*

**CE-IVD**

**SolGent**  
[www.SolGent.com](http://www.SolGent.com)

# About 2019 Novel Coronavirus (2019-nCoV)


Figure 3: Phylogenetic analysis of full-length genomes of 2019-nCoV and representative viruses of the genus Betacoronavirus  
2019-nCoV=2019 novel coronavirus. MERS-CoV=Middle East respiratory syndrome coronavirus. SARS-CoV=severe acute respiratory syndrome coronavirus.


Figure 5: Phylogenetic analysis and homology modelling of the receptor-binding domain of the 2019-nCoV, SARS-CoV, and MERS-CoV

2019-nCoV was closely related with **88% identity** to **two bat-derived SARS-like coronaviruses**, bat-SL-CoVZC45 and bat-SL-CoVZXC21. 2019-nCoV was closely related with **SARS-CoV (about 79%)** and **MERS-CoV (about 50%)**.

Ref: Genomic characterization and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. Published Online January 29, 2020 [https://doi.org/10.1016/S0140-6736\(20\)30251-8](https://doi.org/10.1016/S0140-6736(20)30251-8)


2019 Novel Coronavirus (2019-nCoV) is a virus (more specifically, a [coronavirus](#)) identified as the cause of an outbreak of respiratory illness first detected in Wuhan, China. Early on, many of the patients in the outbreak in Wuhan, China reportedly had some link to a large seafood and animal market, suggesting animal-to-person spread. However, a growing number of patients reportedly have not had exposure to animal markets, indicating person-to-person spread is occurring. At this time, it's unclear how easily or sustainably this virus is spreading between people. The latest situation summary updates are available on CDC's web page [2019 Novel Coronavirus, Wuhan, China](#). (Ref: <https://www.cdc.gov/coronavirus/2019-ncov/about/index.html>)

# WHO COVID-19(2019-nCoV) technical guidance

## Laboratory testing for 2019-nCoV in humans

Country	Institute	Gene targets
China	China CDC	ORF1ab and N
Germany	Charité	RdRP, E
Hong Kong	HKU	ORF1b-nsp14, N
Japan	National Institute of Infectious Diseases	ORF1a
Thailand	National Institute of Health	N
US	US CDC	Three N primers, Rnase P
KOREA	KCDC	RdRp, E

# Comparison of COVID-19 Detection region by country

**SolGent**

2019 Novel Coronavirus  
Detection system



*DiaPlexQ™* 2019 Novel Coronavirus  
(2019-nCoV) Detection kit

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
SolGent Detection target region 1

N gene

† 2019-Novel Coronavirus (2019-nCoV) Real-time qRT-PCR Panel Primers and Probes by CDC

\* National Institute for Viral Disease Control and Prevention

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

## ► Instrument

- Applied Biosystems™ 7500, 7500 Fast Instrument System
- Bio-Rad CFX96™ System

# DiaPlexQ™ *Novel Coronavirus (2019-nCoV) Detection Kit*


## Amplification Plot information

## Fluorescence information

Target	5' Fluorophore	3' Quencher
<i>N gene</i>	FAM	BHQ1
<i>ORF1a</i>	VIC / JOE*	BHQ1
PCR control	Texas Red/ Cal Fluor Red 610*	BHQ2

\*ABI 7500 / 7500 Fast: JOE, Texas Red | Bio Rad CFX96™: VIC, Cal Fluor Red 610

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## Sample collection

Sample collection according to sample type

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## Multiplex OneStep RT-qPCR

*DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit*

## Data analysis

Use of software for each instrument

- HotStart PCR: high-specificity
- OneStep PCR : multiple targets in a single reaction
- Reliable system: automatic PCR control
- Easy-to-use master mix: just adding template and Primer/Probe Mix
- Rapid detection: OneStep RT-qPCR

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## Sample collection

Sample collection according to sample type

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## Multiplex OneStep RT-qPCR

**DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit**

## Data analysis

Use of software for each instrument

Cap color	Reaction Mixture	Vol.
Red ●	2X OneStep RT-qPCR Buffer	10 $\mu$ l
Blue ●	OneStep RT-qPCR Enzyme mix	2 $\mu$ l
Violet ●	Primer/Probe Mixture	3 $\mu$ l
	Sample template	5 $\mu$ l
	<b>Total</b>	<b>20 <math>\mu</math>l</b>

PCR Condition	Reverse Transcription	50 °C	15 min	X 1
	Initial PCR activation	95 °C	15 min	X 1
	Denaturation	95 °C	20 sec	X 45
	Annealing / Extension	60 °C	40 sec	

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- PC should be **positive** and with CT values within 35 cycles. If PC are *negative*, the testing results for that plate are invalid. Repeat qRT-PCR test.
- NTC should be **negative**. If NTCs are *positive*, the testing results for that plate are invalid.
- Detection of the PCR Control(PCRC) in the Texas Red(or Cal Fluor Red 610) detection channel is not required for positive results in the FAM or JOE(VIC) detection channel. A high Virus RNA load in the sample can lead to reduced or absent PCR Control signal.
- If the results show  $40 < Ct \leq 45$ , perform the experiment again.

\* RT-PCR inhibition or reagent failure. Repeat testing from original sample or collect and test a new sample.

# Certifications of DiaPlexQ™ Novel Corona Detection kit

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Director of High-Tech Med. Dept. of Medical Dev. National Institute of Health Evaluation, Ministry of Food and Drug Safety

Osong Health Technology, Ad 167 Osongsaengmyeong2-ro, Osong (Seongju-si), Chungcheong Province, Korea Tel: +82-43-719-2396

No. of Certificate : 2020022773

### Certificate of Free Sales

Exporting(certifying) country : Republic of Korea  
Importing(requesting) country : Thailand

The Ministry of Food and Drug Safety, certifies that the following fire is manufacture medical devices under the Medical Device Act and the following permitted to be freely sold in overseas markets.

Manufacturer (Registered No. : 4010) : SolGent Co., Ltd.

1F, 2F, 43-10, Techno 5-ro, Yuseong-gu, Daejeon

Product-Licence No.	Classification
20-140	ISO 9001 to ISO 13485 (Quality Management System) and ISO 14971 (Risk Management) for Medical Devices

Product Name : DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

Registered List of Product Classification and Code

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2. The medical device(s) is/are not covered by Annex II and self-testing according to Directive 98/79/EC, with the applicable requirements of the following documents

Ref. No.  
EN ISO 15223-1:2016  
ISO 13485:2016  
ISO 13612:2002  
EN ISO 17511:2003  
EN 23840:2015  
EN 13641:2002  
EN ISO 14971:2012  
EN ISO 18113-1:2011  
EN ISO 18113-2:2011  
IEC 62366-1:2015  
IEC 62366-2:2016

3. The medical device(s) is/are not covered by Annex II and self-testing according to Directive 98/79/EC, with the applicable requirements of the following documents

Ref. No.  
EN ISO 15223-1:2016  
ISO 13485:2016  
ISO 13612:2002  
EN ISO 17511:2003  
EN 23840:2015  
EN 13641:2002  
EN ISO 14971:2012  
EN ISO 18113-1:2011  
EN ISO 18113-2:2011  
IEC 62366-1:2015  
IEC 62366-2:2016

4. The medical device(s) is/are not covered by Annex II and self-testing according to Directive 98/79/EC, with the applicable requirements of the following documents

Ref. No.  
EN ISO 15223-1:2016  
ISO 13485:2016  
ISO 13612:2002  
EN ISO 17511:2003  
EN 23840:2015  
EN 13641:2002  
EN ISO 14971:2012  
EN ISO 18113-1:2011  
EN ISO 18113-2:2011  
IEC 62366-1:2015  
IEC 62366-2:2016

February, 2020 Signature: *Chamajong*

SolGent Co., Ltd. Solutions for Genetic Technologies  
SolGent Co., Ltd. 02/2020 V2.0


CFS of KFDA

Approved at third by KCDC-EUA

CE-IVD

# DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

## The Analytical Sensitivity Test by viral RNA




Viral RNA Concentration	Total Copies / rxn	Hit Rate (%)
63 fg/ul	10,000	100
6.3 fg/ul	1,000	100
630 ag/ul	100	100
63 ag/ul	10	95

\*Viral RNA reference: National Hospital Resource Bank (NCCP No. 43326)

The analytical sensitivity of the DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit determined by Probit analysis is 100 copies/ul of viral RNA

# DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

The Analytical Sensitivity Test by In vitro synthesized **Transcripts RNA**


Transcripts RNA Conc. of N gene	Transcripts RNA Conc. of Orf1a	Total Copies / rxn	Hit Rate (%)
7.5 fg/ul	1.25 fg/ul	10,000	100
750 ag/ul	125 ag/ul	1,000	100
75 ag/ul	12.5 ag/ul	100	100
7.5 ag/ul	1.25 ag/ul	10	95

The analytical sensitivity of the DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit determined by Probit analysis is 10 copies/ul of in vitro transcribed RNA


# DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

## The Analytical Sensitivity Test by Target region Plasmid DNA


N gene – Single Plasmid DNA

ORF1a – Single Plasmid DNA

ORF1a & N gene – Multi Plasmid DNA


Plasmid DNA Conc. of N gene	Plasmid DNA Conc. of Orf1a	Total Copies / rxn	Hit Rate (%)
200 fg/ul	200 fg/ul	10,000	100
20 ag/ul	20 ag/ul	1,000	100
2 ag/ul	2 ag/ul	100	100
0.2 ag/ul	0.2 ag/ul	10	95

The analytical sensitivity of the DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit determined by Probit analysis is 10 copies/ul of in vitro transcribed RNA

# DiaPlexQ™ *Novel Coronavirus (2019-nCoV)* Detection Kit


Organisms tested to demonstrate the analytical *Specificity*  
(*Respiratory virus associated, Human RNA*)

**Control Template (2019-nCoV)**

- N gene (FAM)
- ORF1a (JOE)
- PCRC (Texas Red)

**38 negative strains**

- N gene (FAM)
- ORF1a (JOE)
- PCRC (Texas Red)


Organisms	FAM or JOE (N or Orf1a)	Texas Red (PCRC)
Parainfluenza I	N/D	≤ 26
Parainfluenza II	N/D	≤ 26
Parainfluenza III	N/D	≤ 26
Parainfluenza IV	N/D	≤ 26
Influenza A	N/D	≤ 26
Influenza B	N/D	≤ 26
Adenovirus	N/D	≤ 26
Respiratory syncytial virus A	N/D	≤ 26
Respiratory syncytial virus B	N/D	≤ 26
Rhino B, A	N/D	≤ 26
Bocavirus	N/D	≤ 26
Metapneumovirus	N/D	≤ 26
Beta Coronavirus OC43	N/D	≤ 26
Alpha Coronavirus 229E	N/D	≤ 26
Enterovirus	N/D	≤ 26
Human total RNA (10ng/μl)	N/D	≤ 26
Control Template (2019-nCoV)	N/D	≤ 26

*DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit did **not Cross-react** with any of specified organisms*

# DiaPlexQ™ **Novel Coronavirus (2019-nCoV) Detection Kit**


Organisms tested to demonstrate the analytical **Specificity**  
(Pneumonia and Tuberculosis associated)

**Control Template (2019-nCoV)**

- N gene (FAM)
- ORF1a (JOE)
- PCRC (Texas Red)

**38 negative strains**

- N gene (FAM)
- ORF1a (JOE)
- PCRC (Texas Red)

Organisms	FAM or JOE ( N or Orf1a)	Texas Red (PCRC)
<i>Acinetobacter baumannii</i>	N/D	≤ 26
<i>Bordetella parapertussis</i>	N/D	≤ 26
<i>Bordetella pertussis</i>	N/D	≤ 26
<i>Chlamydomphila pneumoniae</i>	N/D	≤ 26
<i>Haemophilus influenza</i>	N/D	≤ 26
<i>Klebsiella pneumoniae</i>	N/D	≤ 26
<i>Legionella pneumophila</i>	N/D	≤ 26
<i>Moraxella catarrhalis</i>	N/D	≤ 26
<i>Mycoplasma pnemoniae</i>	N/D	≤ 26
<i>Pseudomonas aeruginosa</i>	N/D	≤ 26
<i>Serratia marcescens</i>	N/D	≤ 26
<i>Staphylococcus aureus</i>	N/D	≤ 26
<i>Stenotrophomonas maltophilia</i>	N/D	≤ 26
<i>Streptococcus pneumoniae</i>	N/D	≤ 26
<i>Mycobacterium abscessus</i>	N/D	≤ 26
<i>Mycobacterium avium</i>	N/D	≤ 26
<i>Mycobacterium bovis</i>	N/D	≤ 26
<i>Mycobacterium chelonae</i>	N/D	≤ 26
<i>Mycobacterium intracellulare</i>	N/D	≤ 26
<i>Mycobacterium kansasii</i>	N/D	≤ 26
<i>Mycobacterium scrofulaceum</i>	N/D	≤ 26
<i>Mycobacterium tuberculosis</i>	N/D	≤ 26

**DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit did **not Cross-react** with any of specified organisms**

- Ordering information for relative Products

Name of kits	Cat. No.	Detection Method
<i>DiaPlexQ™</i> Novel Coronavirus (2019-nCoV) Detection Kit	SQD52-K100 SQD52-K020	Real-Time PCR
<i>DiaPlexQ™</i> 2019-nCoV (RdRp, E, N) Detection Kit	SQD55-K100 SQD55-K020	
<i>DiaPlexQ™</i> RV16 Detection Kit	SQD50-K100 SQD50-K020	
<i>DiaPlexQ™</i> PneumoPatho16 Detection Kit	SQD80-K100 SQD80-K020	
<i>DiaPlexQ™</i> MTC/NTM detection Kit-Ver. 3.0	SQD25-K100 SQD25-K020	

## DiaPlexQ™ RV16 Detection Kit

### ▶ Target Organisms

SET/01 | *Parainfluenza virus I*  
*Parainfluenza virus II*  
*Parainfluenza virus III*

SET/02 | *Parainfluenza virus IV*  
*Influenza A virus*  
*Influenza B virus*

SET/03 | *Adenovirus*  
*Respiratory syncytial virus A/B*  
*Rhinovirus*

SET/04 | *Enterovirus*  
*Bocavirus*  
*Metapneumovirus*

SET/05 | *Alpha Corona virus 229E/NL63*  
*Beta Corona virus OC43*  
*MERS Corona virus*

## DiaPlexQ™ PP16 Detection Kit

### ▶ Target Organisms

SET/01 | *Mycoplasma pneumoniae (MP)*  
*Klebsiella pneumoniae (KP)*  
*Legionella pneumophila (LP)*


SET/02 | *Streptococcus pneumoniae (SP)*  
*Staphylococcus aureus (SA)*  
*Chlamydomphila pneumoniae (CP)*

SET/03 | *Pseudomonas aeruginosa (PA)*  
*Moraxella catarrhalis (MC)*  
*Bordetella pertussis (BP)*

SET/04 | *Haemophilus influenza (HI)*  
*Acinetobacter baumannii (AB)*  
*M. tuberculosis (TB) / M. avium (AV)*

SET/05 | *Serratia marcescens (SMar)*  
*Bordetella parapertussis (BPara)*  
*Stenotrophomonas maltophilia (SMalt)*

**NEW**

긴급사용승인(KCDC EUAL), CE-IVD

# DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

(Cat.no. SQD52-K100, SQD52-K020)

Multiplex OneStep qRT-PCR based assay system  
Qualitative Detection of 2019-nCoV  
Control Template (2019-nCoV) included

## Detection Target Region

- ORF1a
- N gene

## Specimen type

- Nasopharyngeal swab
- Oropharyngeal swab
- Sputum

## Instruments

- Applied Biosystems™ 7500 Real-Time PCR Instrument System
- Applied Biosystems™ 7500 Fast Real-Time PCR Instrument System
- Bio-Rad CFX96™ System

## Process

### Sample collection

- Sample collection according to sample type

### Nucleic acid isolation

- Manual Method
- Auto extraction Method


### Multiplex qPCR

DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit

### Data analysis

- Use of software For each instrument

## Result Analysis

Amplification Plot information

Target	5' Fluorophore	3' Quencher
N gene	FAM	BHQ1
ORF1a	VIC / JOE*	BHQ1
PCR control	*Texas Red/ Cal Fluor Red 610	BHQ2

\*ABI 7500 / 7500 Fast: JOE, Texas Red  
Bio Rad CFX96™ : VIC, Cal Fluor Red 610

## Ordering Information for Relative Products

Name of Kits	Cat. No.	Approval	Detection Method
DiaPlexQ™ Novel Coronavirus (2019-nCoV) Detection Kit	SQD52-K100	KCDC EUAL, CE-IVD	Real-Time PCR
DiaPlexQ™ 2019-nCoV (RdRp, E, N) Detection Kit	SQD55-K100	CE-IVD	
DiaPlexQ™ RV16 Detection Kit	SQD50-K100	CE-IVD	
DiaPlexQ™ PneumoPatho 16 Detection Kit	SQD80-K100	RUO	
DiaPlexQ™ MTC/NTM Detection Kit - Ver 3.0	SQD25-K100	CE-IVD	

# SolGent Co.,Ltd.

## DiaPlexQ™ Novel Coronavirus Detection Kit


### Features

1. Multiplex OneStep qRT-PCR
2. Hot Start PCR system by using optimized Hot Start polymerase
3. Commercial Real-time PCR Instrument available
4. High specificity: simultaneous detection of ORF1a and N gene.

	Components	SQD52-K020	SQD52-K100
	2X OneStep qRT-PCR Buffer (2019-nCoV)	200 $\mu$ l x 1 ea	1.0 mL x 1ea
OneStep qRT-PCR Enzyme mix (2019-nCoV)	40 $\mu$ l x 1 ea	200 $\mu$ l x 1 ea	
Primer & Probe Mixture (2019-nCoV)	60 $\mu$ l x 1 ea	300 $\mu$ l x 1 ea	
Control Template (2019-nCoV)	20 $\mu$ l x 1 ea	100 $\mu$ l x 1 ea	

### Detection Targets

#### ■ Simultaneous Detection of ORF1a / N gene

- CDC 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes
- High specific targets were selected based on the Chinese CDC and US CDC.

## Performance

- Analytical sensitivity is ~10 copies. Sensitivity is measured using the RNA(In vitro transcript) with the target gene region of Novel Coronavirus (2019-nCoV) as a template.
- It showed high specificity and not detected on other viruses and bacteria other than new coronavirus as a result of specificity test for respiratory virus, pneumonia, and 38 species of tuberculosis.

## Related Products (16 species of respiratory virus detection kit including like-coronavirus)

1. The Kit consists of 5 SETs, Including MERS-CoV, which detects various respiratory viruses.
2. Coronavirus OC43/229E/NL63 vs. MERS-CoV(simultaneous detection of upE/ORF1a)

SET I	SET II	SET III	SET IV	SET V
<i>Parainfluenza virus I</i>	<i>Influenza virus A</i>	<i>Adenovirus</i>	<i>Enterovirus</i>	<i>Coronavirus OC43</i>
<i>Parainfluenza virus II</i>	<i>Influenza virus B</i>	<i>Respiratory Syncytial virus A/B</i>	<i>Bocavirus</i>	<i>Coronavirus 229E/NL63</i>
<i>Parainfluenza virus III</i>	<i>Parainfluenza virus IV</i>	<i>Rhinovirus A/B/C</i>	<i>Metapneumovirus</i>	<i>Coronavirus MERS</i>

Product Marketing Manager

Business support team: Junho Park (C.P : 010-7541-1525, Tel : 070-7893-7847, E.mail junho.park@solgent.com)

**SolGent Co.,Ltd.**

Tel. +82-1544-5695 | Address: 43-10, Techno 5ro, Yuseong-gu, Daejeon, 34014, Korea

# EC DECLARATION OF CONFORMITY

## SolGent Co., Ltd.

Head Office : 3F, 32, Techno 6-ro, Yuseong-gu, Daejeon, 34014, Korea  
 Factory : 1F, 2F, 43-10, Techno 5-ro, Yuseong-gu, Daejeon, 34014, Korea  
 Tel : +82-42-864-5695, Fax: +82-42-864-5690, global@solgent.com

Declares that the medical device(s) described hereafter  
Other Virology – NA Reagents, 15 04 40 90 00 (EDMA code)

Model Name: DiaPlexQ™ 2019-nCoV (RdRp, E, N) Detection Kit

Catalogue Number: SQD55-K100, SQD55-K020

Has been classified as Others not covered by Annex II and self-testing according to Directive 98/79/EC.

Is in conformity with the applicable requirements of the following documents

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Is subject to the conformity assessment procedure set out in Annex III of Directive 98/79/EC

26th, February, 2020 Signature: Chameejeong

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### Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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**DiaPlexC RV13 Detection Kit is designed to detect Metapneumovirus (MPV), Coronavirus OC43(CoV OC43), Coronavirus 229E (CoV 229E), Enterovirus (EntV), Parainfluenza I, II, III (PIV-I, II, III), Respiratory syncytial virus (RSV), Adenovirus (AdV), Bocavirus (Bov), Rhinovirus A/B/C (RV A/B/C) and Influenza virus A/B (Inf A/B) using multiplex PCR technology. DiaPlexQ Ebola Virus Detection Kit is used for the qualitative detection of Ebola Virus from blood or serum of Ebola Virus infected patients by OneStep qRT-PCR. DiaPlexQ RV16 Detection Kit ualitative detection of Parainfluenza-I, Parainfluenza-II, Parainfluenza-III, Parainfluenza-IV, Influenza A (Inf A), Influenza B (Inf B), Adenovirus (AdV), Respiratory syncytial virus (RSV A & B), Rhinovirus (RV A/B/C), Enterovirus (EntV), Bocavirus (BoV), Metapneumovirus (MPV), Beta Coronavirus OC43 (CoV OC43), Alpha Coronavirus 229E (CoV 229E), Alpha Coronavirus NL63 (CoV NL63), MERS Coronavirus (MERS)) from sputum or Nasopharyngeal smear samples of each Virus infected patients by OneStep qRT-PCR. LytePrime Mplex ZCD RT-PCR Kit 1.0 (ZIKV/CHIKV/DENV) is intended for qualitative detection of Zika virus (ZIKV), Chikungunya virus (CHIKV) and Dengue virus (DENV) from human serum and whole blood. DiaPlexQ ZCD Detection Kit (ZIKV, CHIKV, DENV) is designed to detect Zika virus, Chikungunya virus and Dengue virus rapidly and accurately. DiaPlexQ Novel Coronavirus (2019-nCoV) Detection Kit is in vitro diagnostic reagent to qualitative detection of ORF1a, N gene of COVID-19 by Multiplex OneStep qRT-PCR. DiaPlexQ 2019-nCoV (RdRp, E, N) Detection Kit is in vitro diagnostic reagent to qualitative detection of ORF1ab(RdRp), E gene, N gene of COVID-19 by Multiplex OneStep qRT-PCR. For professional use only. Internal Admin.-No.: SOL-08, SOL-08-01 to SOL-08-05.**

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Bearbeiter / Person responsible  
**Frau Caroline Bauer**

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3F, 32, Techno 6-ro, Yuseong-gu, Daejeon, 34014, Rep. of KOREA

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1F, 2F, 43-10, Techno 5-ro, Yuseong-gu, Daejeon, 34014, Rep. of KOREA

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(We hereby certify that the above manufacturer complies with Korea Good Manufacturing Practices of Medical Devices for the product group listed above)

발행일자(Date of Issue) : 2019. 09. 24

유효기간(Date of Expiration) : 2022. 10. 22

대전지방식품의약품안전청장  
DAEJEON REGIONAL FOOD AND DRUG ADMINISTRATION

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## PRECAUTION FOR TREATMENT AND STORAGE IN DELIVERY

2020 March 06<sup>th</sup>

I. This diagnostic kit will be shipped with an ice pack for the product quality on the stability

※ For custom clearance, packing the kit in dry ice for transport is not possible.

II. Multiple freeze and thaw may effect performance of product, please refrain freeze-thaw.

III. The retention period of the product quality is 5 days, if keeps under refrigeration.

IV. The shipping box should be avoided exposure to direct sunlight and heat on its way.


IV. The shipping box should be required refrigeration(-20°C) before the customer receive.

※ We have no obligation regarding any claim based on failure to comply with the above requirement.

Yours Sincerely,

Jae Hyung You & Du-su Seock / CEO of SolGent Co., Ltd.

*Jay you Andrew seock*

Date : March 06<sup>th</sup> , 2020

www.solgent.com

**THANK YOU**

**SolGent** | [www.SolGent.com](http://www.SolGent.com)

Copyright © 2020 SolGent Co., Ltd. All Rights Reserved.  
This PowerPoint file is SolGent's knowledge property. To distribute, copy, publish is not allowed without permission. If you want to use this file, please contact to SolGent.